

## Lyoph-Ready CesiumTaq

Cat #: GF200



**Amount:** 4000 x 25 µl reactions up to 1 kb (equivalent to 200ul standard enzyme. Volume may be up to 2.5x higher)

**Shipping conditions:** Ice Pack

**Storage conditions:** 4°C for 4 months or -20°C for 2 years with up to 10 freeze/thaw cycles

**Thermostability:** Retains at least 85% activity after 1 hour at 95°C

**Expiration:** On tube label

### PRODUCT DESCRIPTION:

A lyoph-ready preparation of CesiumTaq, a double cold-sensitive mutant of Taq DNA polymerase. Due to its suppressed activity at low temperatures this enzyme is designed for hot-start PCR performance. 10x buffer composition is: 500 mM Tris-Cl pH 8.3, 160 mM ammonium sulfate, 0.25% Brij 58, and 25 mM magnesium chloride.

### TYPICAL PCR PROTOCOL for a 25 µl reaction:

| Reagent                              | Volume                       | Final Concentration |
|--------------------------------------|------------------------------|---------------------|
| 10x Taq Mutant Reaction Buffer       | 2.5 µl                       | 1x                  |
| dNTP mix (10 mM each)                | 0.5 µl                       | 200 µM each         |
| Left Primer                          | variable                     | 200 nM              |
| Right Primer                         | variable                     | 200 nM              |
| DNA template†                        | variable                     | 0.1-100 ng          |
| PCR Enhancer Cocktail (recommended)* | 12.5 µl                      | 1x                  |
| CesiumTaq                            | 0.05 – 0.25 µl**             |                     |
| De-ionized distilled H2O             | Adjust final volume to 25 µl |                     |

† DNA amount depends mostly on genome size and target gene copy number.

\* For optimal performance, we recommend using one of our PCR Enhancer Cocktails (PEC-1, PEC-1GC, PEC-2, or PEC-2-GC) or 1.3 M Betaine, a generic PCR enhancer.

\*\* To determine specific optimal enzyme concentration, we strongly recommend an enzyme titration test for each target. A good starting amount of the enzyme per 25 µl reaction is 0.05 µl. Targets larger than 1 kb may require more enzyme or may benefit from the use of an LA (Long Accurate) version of the polymerase.

### CYCLING CONDITIONS:

1. Denaturing: 94° for 2-8 minutes for 1 cycle
2. Denaturing: 94° for 40-60 seconds
3. Annealing: 50°-68° depending on the specific Tm primers for 40-60 seconds
4. Extension: 68° for 2 min/kb target
5. Repeat steps 2-4 for 25-40 cycles

**Please visit us on the web at [www.klentaq.com](http://www.klentaq.com) for troubleshooting and detailed protocols.**

### REFERENCES:

Kermekchiev, M.B., et al. (2003) Cold-sensitive mutants of Taq DNA polymerase provide a hot start for PCR. Nucl Acids Res. 31, 6139-6147.